

## Roadster S3 Specifications

---

Product Number	S731RS3
Weight Capacity	300 lbs
Drivetrain	RWD
Seat Width	18"
Seat Depth	18"
Backrest Height	14.5"
Basket Size	12.5"W x 10"H x 6"D
Seat-to-Ground Height	19" - 21"
Seat-to-Deck Height	16" - 18"
Overall Length	40.5"
Overall Width	25.5"
Overall Height	36.6"
Turning Radius	47"
Ground clearance	2"
Weight of Heaviest Piece	29 lbs
Total Weight	103 lbs
Max Speed	4 mph
Per Charge Range	up to 12 miles
Motor	DC 24V, 180W
Controller	PG S-DRIVE 45A
Battery Size (2)	12V / 12AH
Charger	1.8A off-board
Gradient	6°
Front Wheel(s)	7" PU
Rear Wheels	8" PU
HCPCS Code	n/a
Standard Colors:	Cherry Red or Galaxy Blue

