

merits

# ROADSTER S3

S731RS3

Three Wheel Scooter w/Lift Option



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features

- Tiller angle adjustable for ergonomic driving position and easy on and off.
- Loop tiller with throttle lever to maneuver by palm or finger.
- Transport and store easily with auto latching lockup design.
- Detachable battery pack with built in charging port.
- Large plastic basket.
- Easily disassembles into 5 pieces.
- Extra charging port under dashboard for easy plug in.
- Simple dashboard easy to read.
- No tools needed for set-up.
- High ground clearance of 2"
- Built-in charging port
- Headlight

## Options

- 40" seat belt
- Multi-Accessory Receiver Adapter
- Oxygen tank holder
- Oxygen tank bag holder
- Crutch holder
- Cane holder
- Walker holder
- (Lg) Foldable rear basket
- (Med) Foldable rear basket
- Cup holder
- Safety flag
- Cover

## Specifications

Product#	S731RS3
Weight Capacity	300 lbs
Drivetrain	RWD
Maximum Speed	4 mph
Per-Charge Range	up to 12 mi.
Battery Size (2 required)	12V/12AH
Ground Clearance	2"
Turning Radius	47"
Overall Length	40.5"
Overall Width	25.5"
Overall Height	36.6"
Seat-to-Ground Height	19" - 21"
Seat-to-Deck Height	16" - 18"
Seat Width	18"
Seat Depth	18"
Backrest Height from seat	14.5"
Basket Size	12.5"W x 10"H x 6"D
Total Weight	103 lbs.
Motor	24V 180W
Controller	PG S-DRIVE 45A
Charger	1.8A off-board
Gradient	6°
Front Wheels	7" PU
Drive Wheels	8" PU
HCPCS Code	n/a
Color Options	<div style="display: flex; align-items: center;"> <div style="width: 15px; height: 15px; background-color: red; margin-right: 5px;"></div> Cherry Red         </div> <div style="display: flex; align-items: center; margin-top: 5px;"> <div style="width: 15px; height: 15px; background-color: blue; margin-right: 5px;"></div> Galaxy Blue         </div>



**merits** *Enhancing life through mobility.*

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com

