

Roadster S4 Specifications

Product Number	S741RS4
Weight Capacity	260 lbs
Drivetrain	RWD
Seat Width	18"
Seat Depth	17"
Backrest Height From Seat	15"
Basket Size	12.5"W x 10"H x 6"D
Seat-to-Ground Height	19" - 21"
Seat-to-Deck Height	16" - 18"
Overall Length	30"
Overall Width	21"
Overall Height	36"
Turning Radius	43.3"
Ground Clearance	2"
Weight of Heaviest Piece	36 lbs
Total Weight	107 lbs
Max Speed	4 mph
Per Charge Range	up to 6 miles
Motor	DC 24V, 100W
Controller	PG S-Drive 45A
Battery Size (2)	12V/12AH
Charger	1.8A off-board
Gradient	6°
Front Wheel(s)	8" PU
Rear Wheels	8" PU
HCPCS Code	n/a
Standard Colors:	Cherry Red and Galaxy Blue

