

merits

ROADSTER S4

S741RS4

Four Wheel Travel-Friendly Scooter



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- Tiller angle adjustable for ergonomic driving position and easy on and off.
- Loop tiller with throttle lever to maneuver by palm or finger.
- Transport and store easily with auto latching lockup design.
- Detachable battery pack with built in charging port.
- Large plastic basket.
- Easily disassembles into 5 pieces.
- Extra charging port under dashboard for easy plug in.
- Simple dashboard easy to read.
- No tools needed for set-up.
- High ground clearance of 2"
- Built-in charging port
- Headlight

Options

- 40" seat belt
- USB charger
- Multi-Accessory Reciever Adapter
- Oxygen tank holder
- Oxygen tank bag holder
- Crutch holder
- Cane holder
- Walker holder
- (Lg) Foldable rear basket
- (Med) Foldable rear basket
- (Lg) Foldable rear basket
- Cup holder
- Safety flag
- Cover

Specifications

Product#	S741RS4
Weight Capacity	260 lbs.
Drivetrain	RWD
Maximum Speed	4 mph
Per-Charge Range	up to 6 mi.
Battery Size (2 required)	12V/12AH
Ground Clearance	2"
Turning Radius	43.3"
Overall Length	30"
Overall Width	21"
Overall Height	36"
Seat-to-Ground Height	19" - 21"
Seat-to-Deck Height	16" - 18"
Seat Width	18"
Seat Depth	17"
Backrest Height from seat	15"
Basket Size	12.5"W x 10"H x 6"D
Total Weight	107 lbs.
Weight of Heaviest Piece	36 lbs.
Motor	DC 24V, 100W
Controller	PG S-DRIVE 45A
Charger	1.8A off-board
Gradient	6°
Front Wheels	8" PU
Drive Wheels	8" PU
HCPCS Code	n/a
Color Options	■ Cherry Red ■ Galaxy Blue



merits *Enhancing life through mobility.*

©Merits Health Products, Inc.
 4245 Evans Avenue Fort Myers, FL 33901
 www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0225

merits