

Silverado Extreme Specifications

Product Number	S941L
Weight Capacity	450 lbs
Drivetrain	RWD
Seat Width	20"
Seat Depth	18"
Backrest Height	25"
Seat-to-Ground Height	26.4"
Seat-to-Deck Height	19"
Overall Length	61"
Overall Width	29"
Overall Height	57.3"
Turning Radius	82.7"
Ground Clearance	5.1"
Weight of Heaviest Piece	185 lbs
Total Weight	387 lbs
Max Speed	9.5 mph
Per Charge Range	up to 55 miles
Motor	DC 24V, 650W
Controller	P&G S-Drive 120A
Battery Size (2)	12V/100AH
Charger	8A off-board
Gradient	10°
Wheel Base	45"
Front Wheel(s)	14" pneumatic
Rear Wheels	16" pneumatic
Brake	Intelligent, Regenerative, Eletromagnetic Brakes
Standard Colors	Black

