

Travel Ease Specifications

Product Number	P101, P181
Weight Capacity	300 lbs, 450 lbs
Drivetrain	RWD
Seat Type	Sling
Seat Width	16" - 22"
Seat Depth	16" - 18"
Backrest Height	15" - 19"
Seat-to-Ground Height	18" - 20"
Overall Length	43" - 45"
Overall Width	22" - 28"
Overall Height	36"
Turning Radius	23" - 25"
Ground clearance	4.8"
Total Weight	146 lbs - 150 lbs
Weight of Heaviest Piece	86 lbs - 97 lbs
Max Speed	4 mph
Battery Size (2)	12V / U1 (34AH); 12V / 22NF (50AH)
Per Charge Range	up to 18 miles
Motor	DC 24V / 180W, DC 24V / 300W
Controller	VSI, VR2 70A
Charger	5A off-board
Gradient	10°
Front Wheels	8" foam-filled
Drive Wheels	12.5" foam-filled
Rear Wheels	1.5" anti-tippers
HCPCS Code	K0822 & K0824
Standard Colors	Silvertone

