

merits

TRAVEL EASE P101

Portable Folding Power Chair

K0822 / K0824



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- Optional weight capacity of 300 or 400 lbs.
- 12.5" rubber drive tires for superior traction
- 4.8" ground clearance for easy obstacle clearance
- Foldable for improved storage and transportation

Accessories

- USB charger
- Oxygen tank holder
- Swing-away joystick
- Cup holder
- Side bag
- Weather cover

Options

- Multiple seat sizes available
- Height adjustable flip-up armrests and swing-away footrests
- Elevating leg rests (E0990)
- 4 Seat cushions for maximum comfort



Specifications

Product#	P101, P181
Weight Capacity	300 lbs., 450 lbs.
Drive train	RWD
Seat Type	Sling
Seat Width	16" - 22"
Seat Depth	16" - 18"
Backrest Height	15" - 19"
Seat-to-Floor Height	18" - 20"
Overall Length	43" - 45"
Overall Width	22" - 28"
Overall Height	36"
Turning Radius	23" - 25"
Ground Clearance	4.8"
Total Weight	146 lbs., 150 lbs.
Wt of Heaviest Pc.	86 lbs. - 97 lbs.
Max Speed	4 mph
Range Per Charge	up to 18 mi.
Motor	DC 24V / 180W, DC 24V / 300W
Controller	VSI, VR2 70A
Battery Size (2 required)	12V / U1 (34AH); 12V / 22NF (50AH)
Charger	5A off-board
Gradient	10°
Front Wheels	8" foam-filled
Drive Wheels	12.5" foam-filled
Rear Wheels	1.5" anti-tippers
HCPCS Code	K0822 K0824
Standard Colors	■ Silvertone



merits *Enhancing life through mobility.*

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com

