

Ultra 3 Actuator Specifications

Product Number	P325
Weight Capacity	300 lbs
Drivetrain	MWD
Seat Width	16-22", adjusted in 1" increments
Seat Depth	16-22", adjusted in 1" increments
Backrest Height	18-28"
Seat-to-Ground Height	Power tilt or pan: 16.5" - 19.5", Captain: 20"-23"
Seat- Elevation (Optional)	12"
Tilt (Optional)	0°-50°
Recline (Optional)	90°-110° fixed
Articulating Footplate	Included
Overall Length	36"
Overall Width	24"
Suspension	Front and Rear
Overall Height	40"
Turning Radius	26"
Ground clearance	2.5"
Total Weight	315-384 lbs
Max Speed	5 mph
Battery Size (2)	34Ah, U1 - 55Ah, 22NF - 55 Ah, 22NF Gel
Per Charge Range	up to 19 miles
Motor	2-Pole/240W
Controller	P&G VR2-A2 60A
Charger	6 Amp
Gradient	6°
Front Wheels	Anti-Tipper
Drive Wheels	14" x 3" foam-filled
Rear Wheels	6" x 2" foam-filled
HCPCS Code	K0835
Standard Colors	Red, Blue, White, and Charcoal Gray

