

## Ultra CST Specifications

---

Product Number	P325-2s C
Weight Capacity	300 lbs
Drivetrain	MWD
Seat Width	16"-22", adjusted in 1" increments
Seat Depth	16"-22", adjusted in 1" increments
Backrest Height	18-28"
Seat-to-Ground Height	Power tilt or pan: 16.5"-19.5", Captain: 20"-23"
Tilt (Optional)	0° - 50°
Recline	90° - 105° (manual)
Overall Length	36"
Overall Width	24"
Overall Height	40"
Suspension	Front and rear
Turning Radius	26"
Ground clearance	2.5"
Total Weight	315-384 lbs
Max Speed	5 mph
Battery Size (2)	35Ah, U1 - 55Ah, 22NF - 55Ah, 22NF GEL
Per Charge Range	up to 19 miles
Motor	2-Pole/240W
Controller	P&G VR2-A2 60A
Charger	6 Amp
Gradient	6°
Front Wheels	Anti-Tipper
Drive Wheels	14" x 3" foam-filled
Rear Wheels	6" x 2" foam-filled
HCPCS Code	K0836
Standard Colors	White, Red, Blue and Charcoal Gray

