

Vector HD Specifications

Product Number	P323
Weight Capacity	450 lbs
Drivetrain	MWD
Seat Width	20-24"
Seat Depth	16-22"
Backrest Height	18-28"
Seat-to-Ground Height	18.5" - 19.5"
Seat- Elevation (Optional)	12"
Overall Length	47"
Overall Width	25"
Overall Height	42"
Turning Radius	27"
Ground clearance	3"
Total Weight	315-416 lbs
Max Speed	5 mph
Battery Size (2)	55 Ah, 22NF Gel - 55Ah, 22NF, AGM
Per Charge Range	up to 12.5 miles
Motor	4-Pole / 200W
Controller	P&G VR2 or R-NET
Charger	6 Amp
Gradient	6°
Front Wheels	8" x 3" foam-filled
Drive Wheels	14" x 3" foam-filled
Rear Wheels	8" x 3" foam-filled
HCPCS Code	K0850 / K0858 / K0862
Standard Colors	Candy Apple Red, Metallic Royal Blue, Platinum Silver, and Metallic Black

