

Velocity Specifications

Product Number	P325
Weight Capacity	300 lbs
Drivetrain	MWD
Seat Width	16-22"
Seat Depth	16-22"
Backrest Height	23"
Seat-to-Ground Height	16.5" - 19.5"
Seat- Elevation (Optional)	12"
Overall Length	36"
Overall Width	24"
Overall Height	40"
Turning Radius	26"
Ground clearance	2.5"
Total Weight	320-404 lbs
Max Speed	6 mph
Battery Size (2)	55 Ah, 22NF Gel - 55Ah, 22 NF AGM
Per Charge Range	up to 18 miles
Motor	2-Pole/275W
Controller	P&G VR2 or R-NET
Charger	6 Amp
Gradient	6°
Front Wheels	6" x 2" - solid
Drive Wheels	14" x 3" - foam-filled
Rear Wheels	6" x 2" - solid
HCPCS Code	K0848 / K0849 / K0856 / K0861
Standard Colors	Green, Metallic Black, Pink, Metallic Royal Blue, Candy Apple Red and Platinum Silver

