



# VELOCITY P325

MWD Powerbase

K0848 / K0849 / K0856 / K0861



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features

- Center of gravity tilt 0° - 50° posterior
- Bright LED lights, included
- Shear reduction recline 90° - 150°
- Tension-adjustable front and rear suspension
- Adjustable seat width and depth from 16" - 22"; depth adjustable in 1" increments
- 24" base width and seat-to-floor height as low as 16.5"
- Rated up to 18 miles per charge

## Accessories

Based on chair configuration)

- USB charger
- Cup holder
- Oxygen tank holder
- Cane & crutch holder
- Foldable rear basket

## Options

- 12" telescoping power seat lift
- Joystick knob choices
- Trunk support pad
- Directional swingaway
- Captain, Pan or AVID powered seat
- Multiple footplate options

## Specifications

Product#	P325
Weight Capacity	300 lbs.
Drivetrain	MWD
Seat Width	16" - 22"
Seat Depth	16" - 22"
Backrest Height	23"
Seat-to-Floor Height	16.5" - 19.5"
Seat Elevation (Optional)	12"
Overall Length	36"
Overall Width	24"
Overall Height	40"
Turning Radius	26"
Ground Clearance	2.5"
Total Weight	320 - 404 lbs.
Max Speed	6 mph
Battery Size (2 required)	55Ah, 22NF Gel, 55Ah, 22NF AGM
Range Per Charge	up to 18 mi.
Motor	2-Pole / 275W
Controller	P&G VR2 or R-NET
Charger	6 Amp
Gradient	6°
Front Wheels	6" X 2" - solid
Drive Wheels	14" X 3" foam-filled
Rear Wheels	6" X 2" - solid
HPCSCCode	K0848 / K0849 K0856 / K0861

Color Options

- Green
- Metallic Black
- Pink
- Metallic Royal Blue
- Candy Apple Red
- Platinum Silver



Optional MK batteries



**merits** *Enhancing life through mobility.*

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0125