

Vision CF Specifications

Product Number	P322
Weight Capacity	300 lbs
Drivetrain	FWD
Seat Type	Captain or Pan
Seat Width	16" - 20"
Seat Depth	16" - 20"
Backrest Height	18" - 20"
Seat-to-Floor Height	18" - 22 "
Overall Length	33"
Overall Width	25"
Overall Height	42"
Turning Radius	22.5"
Ground clearance	2.4"
Total Weight	154 lbs
Weight of Heaviest Piece	62.7 lbs
Max Speed	4 mph
Battery Size (2)	12V / U1 (34AH)
Per Charge Range	up to 18 miles
Motor	DC 24V / 130W
Controller	Dynamic LiNX 40A
Charger	4A off-board
Gradient	6°
Front Wheels	3" anti-tippers
Drive Wheels	9" foam-filled
Rear Wheels	6" foam-filled
HCPCS Code	K0815 / K0816 / K0820 / K0821 / K0822 / K0823
Standard Colors	Candy Apple Red

