

merits

# VISION CF P322

FWD Powerbase

K0815 / K0816 / K0820 / K0821 / K0822 / K0823



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

# Features

- A true FWD with a 22.5" turn radius
- 9" knobby drive tires for superior traction
- Dynamic LiNX controller for convenient Bluetooth programming and diagnostics
- Premium, semi-reclining ergonomic seat up to 135°
- 16" - 20" customizable seat depth in 1" increments

# Accessories

- Unoccupied transit tie down kit
- Cup holder
- USB charger
- Oxygen tank holder
- Walker holder
- Foldable rear basket
- Safety flag
- Side under, arm bag

# Options

- Optional, 16" - 20" seat width
- Elevating leg rests (K0195)
- Captain or Pan powered seat option

# Specifications

Product#	P322
Weight Capacity	300 lbs.
Drive train	FWD
Seat Type	Captain or Pan
Seat Width	16" - 20"
Seat Depth	16" - 20"
Backrest Height	18" - 20"
Seat-to-Floor Height	18" - 22"
Overall Length	33"
Overall Width	25"
Overall Height	42"
Turning Radius	22.5"
Ground Clearance	2.4"
Total Weight	154 lbs.
Wt of Heaviest Pc.	62.7 lbs.
Max Speed	4 mph
Range Per Charge	up to 18 mi.
Motor	DC24V / 130W
Controller	Dynamic LiNX 40A
Battery Size (2 required)	12V / U1 (34AH)
Charger	4A off-board
Gradient	6°
Front Wheels	3" anti-tippers
Drive Wheels	9" foam-filled
Rear Wheels	6" foam-filled
HCPCS Code	K0815 / K0816 K0820 / K0821 K0822 / K0823

Color Option ■ Candy Apple Red



merits



**merits** *Enhancing life through mobility.*

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0225