

Vision Sport Specifications

Product Number	P326E
Weight Capacity	300 lbs
Drivetrain	MWD
Seat Type	Captain or Pan
Seat Width	18" - 20"
Seat Depth	16 - 20"
Backrest Height	18" - 20"
Seat-to-Floor Height	17" - 21"
Overall Length	34"
Overall Width	24"
Overall Height	48.5"
Turning Radius	20.6"
Ground clearance	2"
Total Weight	184 lbs
Weight of Heaviest Piece	90 lbs
Max Speed	5 mph
Battery Size (2)	12V / U1 (34AH)
Per Charge Range	up to 18 miles
Motor	DC 24V / 160W
Controller	Dynamic LiNX 40A
Charger	4A off-board
Gradient	6°
Front Wheels	6" - solid
Drive Wheels	10" - solid
Rear Wheels	6" - solid
HCPCS Code	K0822 / K0823
Standard Colors	Red and Black

