

## Vision Sport with Lift Specifications

---

Product Number	P326D
Weight Capacity	300 lbs
Drivetrain	MWD
Seat Type	Captain or Pan
Seat Width	18" - 20"
Seat Depth	16" - 20"
Backrest Height	20"
Seat-to-Floor Height	20.5" - 29.5"
Seat Elevation	7.5"
Overall Length	34"
Overall Width	24"
Overall Height	48.5"
Turning Radius	20"
Ground clearance	2"
Total Weight	214 lbs
Weight of Heaviest Piece	120 lbs
Max Speed	5 mph
Battery Size (2)	12V / U1 (34AH)
Per Charge Range	up to 18 miles
Motor	DC 24V / 160W
Controller	P&G VR2 60A
Charger	4A off-board
Gradient	6°
Front Wheels	6" - solid
Drive Wheels	10" - solid
Rear Wheels	6" - solid
HCPCS Code	K0830 / K0831
Standard Colors	Red or Black

