



VISION ULTRA G2

Captain seat | Pan seat | Single Power

K0822 / K0823 / K0835



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- Center of gravity tilt 0° - 50° posterior
- Bright LED lights, included
- Shear reduction recline 90° - 150°
- Tension-adjustable front and rear suspension
- Adjustable seat width and depth from 16" - 22"; depth adjustable in 1" increments
- 24" base width and seat-to-floor height as low as 16.5"
- Rated up to 18 miles per charge

Options

- 12" elevating seat
- Power articulating foot platform
- Elevating leg rests or swing-away footrests
- Swing-away joystick
- Oxygen tank holder
- USB charger

Models

Vision Ultra - G2, Standard, Captain P325-2S C

Vision Ultra - G2, Standard, Pan P325- 2S S

Vision Ultra - G2, Standard, Single Power P325-2SSP S



Optional MK batteries



Specifications

Product#	P325
Weight Capacity	300 lbs.
Drivetrain	MWD
Maximum Speed	5 mph
Per-Charge Range	up to 19 mi.
Battery Size (2 required)	35Ah, U1; 55Ah, 22NF, 55Ah, 22NF GEL
Ground Clearance	2.5"
Turning Radius	26"
Overall Length	36"
Overall Width	24"
Suspension	Front and Rear
Seat-to-Floor Height	Power tilt or pan: 16.5"-19.5", Captain: 20"-23"
Seat Width	16"-22", adjusted in 1" increments
Seat Depth	16"-22", adjusted in 1" increments
Backrest Height	18"-28"
Total Weight	315 - 384 lbs.
Motor	2-Pole/240W
Controller	P&G VR2-A2 60A
Charger	6 Amp
Gradient	6°
Front Wheels	Anti-Tipper
Drive Wheels	14" x 3" foam-filled
Casters	6" x 2" foam-filled
Color Options	<input type="checkbox"/> White <input checked="" type="checkbox"/> Red <input checked="" type="checkbox"/> Blue <input checked="" type="checkbox"/> Charcoal Gray



merits Enhancing life through mobility.

©Merits Health Products, Inc.
 4245 Evans Avenue Fort Myers, FL 33901
 www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0225