

merits

# ATLANTIS P710

Super Heavy Duty RWD Powerbase

K0826 / K0827



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features

- Heavy-duty power wheelchair with a weight capacity of 600 lbs.
- Four seat posts provide a stable and durable structure for the frame
- Extra-wide seat and extra-high back provide great comfort for bariatric users
- Super high power in-line motors help riders negotiate 10° uphill ramps with ease
- Adjustable seat depth from 20" - 22"; adjustable in 1" increments

## Accessories

- Oxygen tank holder
- Cane & crutch holder
- Walker holder
- Foldable rear basket
- Cup holder
- Side bag
- USB charger
- Cover

## Options

- Captain or Pan power seat option
- 24" or 26" seat width option
- Elevating legrest available (K0195)
- Swing-away joystick
- Flip-up, height-adjustable, detachable armrests



## Specifications

Product#	P710
Weight Capacity	600 lbs.
Drivetrain	RWD
Seat Type	Captain or Pan
Seat Width	24" & 26"
Seat Depth	20" & 22"
Backrest Height	24"
Seat-to-Floor Height	23" - 24"
Overall Length	44"
Overall Width	27"
Overall Height	50"
Turning Radius	31.5"
Ground Clearance	3"
Total Weight	304 lbs.
Wt of Heaviest Pc.	167 lbs.
Max Speed	5 mph
Range Per Charge	up to 32 mi.
Motor	DC 24V / 300W
Controller	P&G VR2 60A
Battery Size (2 required)	75Ah, M24 AGM
Charger	6A off-board
Gradient	10°
Front Wheels	9" pneumatic
Drive Wheels	12" pneumatic
Rear Wheels	3" anti-tippers
HCPCS Codes	K0826 / K0827
Color Options	<span style="color: red;">■</span> Candy Apple Red <span style="color: blue;">■</span> Metallic Blue



merits



**merits** Enhancing life through mobility.

©Merits Health Products, Inc.  
 4245 Evans Avenue Fort Myers, FL 33901  
 www.meritsusa.com | 239-772-0579 | sales@meritsusa.com

