



# AXCEL NARROW P335

FWD Powerbase

K0856 / K0861



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

# Features

- 22.5" wide front wheel drive chair is easy to maneuver in tight spots, narrow hallways and door frames
- Tight 90° turns
- Navigating over uneven terrain is easy with FWD 14" drive tires, maintaining contact.
- Independent hydraulic, high-performance suspension
- Large 9" rear casters to navigate obstacles and rough terrain
- Bright LED lights, included
- Tool-free, rear battery access
- Power recline adjustments, 85° - 160°
- Powered by MK batteries, rated up to 19 miles per charge

# Accessories

(Based on chair configuration)

- USB charger
- Oxygen tank holder
- Cup holder & mount

# Options

- ASL alternative drive
- Therafin positioning products
- Power articulating footplate
- Push handles
- Proform headrest with EZ-Clamp bracket
- 12" seat elevation

# Specifications

Product#	P335
Weight Capacity	300 lbs.
Drivetrain	FWD
Seat Width	14 - 20"
Seat Depth	14 - 20"
Backrest Height	18 - 22"
Seat-to-Floor Height	17.5 - 19.5"
Seat Elevation (Optional)	12"
Overall Length	39.6"
Overall Width	22.5"
Overall Height	41.5"
Turning Radius	26" w/o footplate
Ground Clearance	3"
Total Weight	350-385 lbs.
Max Speed	5.2 mph
Battery Size (2 required)	50Ah, M50 AGM
Range Per Charge	up to 19 mi.
Motor	2-Pole/222W
Controller	R-NET 90A
Charger	6 Amp
Gradient	7.5°
Front Wheels	Anti-Tipper
Drive Wheels	14"x 3" hybrid gel
Rear Wheels	9" hybrid gel
HPCSCCode	K0856/K0861
Color Options	<ul style="list-style-type: none"> <li><span style="display: inline-block; width: 15px; height: 15px; background-color: black; margin-right: 5px;"></span> Black Noir</li> <li><span style="display: inline-block; width: 15px; height: 15px; background-color: gray; margin-right: 5px;"></span> Dove Gray</li> <li><span style="display: inline-block; width: 15px; height: 15px; background-color: blue; margin-right: 5px;"></span> Sapphire Blue</li> <li><span style="display: inline-block; width: 15px; height: 15px; background-color: red; margin-right: 5px;"></span> Ruby Red</li> <li><span style="display: inline-block; width: 15px; height: 15px; background-color: green; margin-right: 5px;"></span> Camo Green</li> </ul>



Features MK batteries



**merits** Enhancing life through mobility.

©Merits Health Products, Inc.  
 4245 Evans Avenue Fort Myers, FL 33901  
 www.meritsusa.com | 239-772-0579 | sales@meritsusa.com

