

merits

EZ-GO P321

Travel-Friendly Power Chair

K0814



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- Ergonomic mid-back seat for additional back support
- Quick connect feature for easy assembly and transportation
- Turning radius of 23" for easy maneuverability in tight spaces
- Weight capacity, 300 lbs.
- Compact and lightweight design, weighs a total of 113 lbs.

Accessories

- Oxygen tank holder
- Cup holder
- Cane & crutch holder
- Mounting bracket for cup holder
- Walker holder
- Side bag
- Foldable rear basket
- USB charger



Specifications

Product#	P321
Weight Capacity	300 lbs.
Drivetrain	RWD
Seat Type	Stadium
Seat Width	18"
Seat Depth	15"
Backrest Height	17"
Seat-to-Floor Height	18" - 20"
Overall Length	27"
Overall Width	23"
Overall Height	43"
Turning Radius	23"
Ground Clearance	1.8"
Total Weight	113 lbs.
Wt of Heaviest Pc.	54 lbs.
Max Speed	4.5 mph
Range Per Charge	up to 12 mi.
Motor	DC 24V / 90W
Controller	P&G NVR2
Battery Size (2 required)	12Ah, ES 12 AGM
Charger	1.5A off-board
Gradient	6°
Front Wheels	7" - solid
Drive Wheels	8" - foam-filled
Rear Wheels	3" anti-tippers
HCPCSCode	K0814
Color Options	■ Red ■ Black

Features MK batteries



merits



merits Enhancing life through mobility.

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0325