

merits

GEMINI P3017, P3011 w/LIFT

Heavy Duty RWD Powerbase

K0824 / K0825 / K0830 / K0831



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- Ergonomic semi-reclining seat that reclines up to 135°
- 14" rubber drive tires for superior traction
- Dynamic LiNX controller for easy programming and diagnostics
- 4" ground clearance for easy obstacle clearance
- Optional 5" power elevating seat.

Accessories

- Oxygen tank holder
- Cane & crutch holder
- Walker holder
- Foldable rear basket
- Cup holder
- Side bag
- Cover
- USB charger

Options

- Captain or Pan powered seat option
- 20" and 22" seat width option
- Elevating leg rest available (K0195)
- 5" power lift option, available
- Swing-away joystick
- XL footplate



Specifications

| | |
|---------------------------|--|
| Product# | P3017, P3011 w/lift |
| Weight Capacity | 450 lbs., 350 lbs. w/lift. |
| Drivetrain | RWD |
| Seat Type | Captain or Pan |
| Seat Width | 20" - 22" |
| Seat Depth | 20" - 22" |
| Backrest Height | 22" - 24" |
| Seat-to-Floor Height | 22.5" - 23.5" |
| Seat Elevation | 5" - P3011 w/lift |
| Overall Length | 41" |
| Overall Width | 24" |
| Overall Height | 45" |
| Turning Radius | 21" |
| Ground Clearance | 4" |
| Total Weight | 242 lbs., 350 lbs. |
| Wt of Heaviest Pc. | 126 lbs., 146 lbs. |
| Max Speed | 5 mph |
| Range Per Charge | up to 20 mi. |
| Motor | DC 24V / 240W |
| Controller | P&G VR2 60A / Dynamic Linx |
| Battery Size (2 required) | 55Ah, 22NF AGM |
| Charger | 5A off-board |
| Gradient | 10° |
| Front Wheels | 9" - foam-filled |
| Drive Wheels | 14" - foam-filled |
| Rear Wheels | 3" anti-tippers |
| HCPCS Codes | K0824 / K0825 K0830 / K0831 |
| Color Options | ■ Candy Apple Red ■ Metallic Silver Blue |

Features MK batteries



merits Enhancing life through mobility.

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0325