

merits

# PIONEER 10 S341

Heavy Duty 4 Wheel Scooter

K0899



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features


- Hydraulic, adjustable tiller
- 22"W x 20"D x 22"H, sliding, captain seat
- Half speed inhibit switch
- Headlight, tail lights, turn signals, and hazard lights
- 4" ground clearance
- Four wheels provide more stability

## Accessories

- Oxygen Tank Holder
- Cane & Crutch Holder
- Walker Holder
- (Lg) Foldable Rear Basket
- (Med) Foldable Rear Basket
- Cover



## Specifications

Product#	S341
Weight Capacity	500 lbs.
Drivetrain	RWD
Maximum Speed	6 mph
Per-Charge Range	up to 32 miles
Battery Size (2 required)	75Ah, M24 AGM
Ground Clearance	4"
Turning Radius	62"
Overall Length	59"
Overall Width	27"
Overall Height	47"
Suspension	Yes
Seat-to-Floor Height	16" - 19"
Seat-to-Ground Height	25" - 28"
Seat Width	22"
Seat Depth	20"
Backrest Height	22"
Total Weight	269 lbs.
Weight of Heaviest Pc.	121 lbs.
Motor	DC 24V, 400W
Controller	Dynamic Rhino2 160A
Charger	6A off-board
Gradient	12°
Front Wheels	12.5" flat-free
Drive Wheels	12.5" flat-free
HCPCS Code	K0899
Color Option	 Red

merits



 **merits** *Enhancing life through mobility.*

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0225