

merits

# PIONEER 3

S131

3 Wheel Scooter

K0807



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features

- Hydraulic, adjustable tiller
- 20"W x 20"D x 20"H, sliding, captain seat
- Half speed inhibit switch
- Headlight, tail lights, turn signals and hazard lights
- 4" ground clearance
- Three wheels provide a tighter turning radius for increased maneuverability

## Accessories

- Oxygen Tank Holder
- Cane & Crutch Holder
- Walker Holder
- (Lg) Foldable Rear Basket
- (Med) Foldable Rear Basket
- Cover



## Specifications

Product#	S131
Weight Capacity	400 lbs.
Drivetrain	RWD
Maximum Speed	5 mph
Per-Charge Range	up to 18 miles
Battery Size (2 required)	33Ah, U1 AGM
Ground Clearance	4"
Turning Radius	39"
Overall Length	47"
Overall Width	23"
Overall Height	45.5"
Suspension	Yes
Seat-to-Deck Height	16" - 19"
Seat-to-Ground Height	25" - 28"
Seat Width	20"
Seat Depth	20"
Backrest Height	20"
Total Weight	178 lbs.
Weight of Heaviest Pc.	86 lbs.
Motor	DC 24V, 250W
Controller	Dynamic Rhino2 90A
Charger	5A off-board
Gradient	8°
Front Wheels	10" flat-free
Drive Wheels	10" flat-free
HCPCS Code	K0807
Color Option	<span style="color: red;">■</span> Red

Features MK batteries



**merits** Enhancing life through mobility.

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



merits

rev. 0225