

merits

# PIONEER 4 S141

Four Wheel Scooter w/Lift Option

K0899



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features

- One-piece frame holds up to 400 lbs.
- Steel bumper comes standard
- Adjustable captain seat
- Padded foam armrests are height, weight, and angle adjustable
- Hazard lights
- Turn signals
- Headlight and tail lights
- Anti-tippers
- Range of up to 18 miles per charge
- 4" ground clearance
- Four wheels provide a tighter turning radius for increased maneuverability
- Optional seat lift

## Accessories

- Oxygen tank holder - attaches to seat bracket
- Cane & crutch holder - attaches to seat bracket
- Walker holder - attaches to seat bracket
- Foldable rear basket (medium or large) - attaches to seat bracket

## Specifications

Product#	S141
Weight Capacity	400 lbs
Drivetrain	RWD
Maximum Speed	5 mph
Per-Charge Range	up to 18 miles
Battery Size (2 required)	33Ah, U1 AGM
Ground Clearance	4"
Turning Radius	52"
Overall Length	48"
Overall Width	23"
Overall Height	45.5"
Seat-to-Ground Height	25" - 28"
Seat-to-Deck Height	16" - 19"
Seat Width	20"
Seat Depth	20"
Backrest Height	20"
Total Weight	189 lbs
Motor	DC 24V, 250W
Controller	Dynamic Rhino2 90A
Charger	5A off-board
Gradient	8°
Front Wheels	10" flat-free
Drive Wheels	10" flat-free
HCPCS Code	K0899
Color Options	<div style="display: inline-block; width: 15px; height: 15px; background-color: red; margin-right: 5px;"></div> Red <div style="display: inline-block; width: 15px; height: 15px; background-color: blue; margin-right: 5px;"></div> Blue



Features MK batteries



merits



**merits** Enhancing life through mobility.

©Merits Health Products, Inc.  
 4245 Evans Avenue Fort Myers, FL 33901  
 www.meritsusa.com | 239-772-0579 | sales@meritsusa.com

