

merits

PIONEER 9 S331

Heavy Duty Three Wheel Scooter

K0808



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- Hydraulic, adjustable tiller
- 22"W x 20"D x 22"H, sliding, captain seat
- Half speed inhibit switch
- Headlight, tail lights, turn signals, and hazard lights
- 4" ground clearance
- Three wheels provide a tighter turning radius for increased maneuverability

Accessories

- USB Charger
- Oxygen Tank Holder
- Oxygen Tank Bag Holder
- Cane & Crutch Holder
- Walker Holder
- (Lg) Foldable Rear Basket
- (Med) Foldable Rear Basket
- Bracket for Cup Holder
- Canopy
- Cover
- Side Bag
- Safety Flag

Specifications

Product#	S331
Weight Capacity	500 lbs
Drivetrain	RWD
Maximum Speed	6 mph
Per-Charge Range	32 miles
Battery Size (2 required)	75Ah, M24 AGM
Ground Clearance	4"
Turning Radius	50"
Overall Length	54"
Overall Width	27"
Overall Height	47"
Suspension	Yes
Seat-to-Floor Height	25" - 28"
Seat-to-Ground Height	25" - 28"
Seat Width	22"
Seat Depth	20"
Backrest Height	22"
Total Weight	249 lbs.
Weight of Heaviest Pc.	101 lbs.
Motor	DC 24V, 400W
Controller	Dynamic Rhino2 160A
Charger	6A off-board
Gradient	12°
Front Wheels	10" flat-free
Drive Wheels	12.5" flat-free
HCPCS Code	K0808
Color Option	■ Red



merits



merits *Enhancing life through mobility.*

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0225