

merits

ROADSTER 3 S731A

Three Wheel Scooter

K0800



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- Tiller angle adjustable for ergonomic driving position and easy transferring on and off the unit
- Loop tiller with throttle lever to maneuver by palm or finger
- Transport and store easily with auto latching lockup design
- Detachable battery pack with built in charging port
- Easily disassembles into 5 pieces
- Extra charging port under dashboard for easy plug in
- No tools needed for set-up
- Adjustable seat height and armrest width
- Rear anti-tip wheels prevent a tip-over backwards

Options

- 40" seat belt
- Multi-accessory receiver adapter
- Oxygen tank holder
- Oxygen tank bag holder
- Crutch and cane holder
- Walker holder
- Medium or Large foldable rear basket
- Cup holder
- Safety flag
- Cover



Specifications

Product#	S731 A
Weight Capacity	300 lbs
Drivetrain	RWD
Maximum Speed	4 mph
Per-Charge Range	up to 6 miles
Battery Size (2 required)	12Ah, ES12 AGM
Ground Clearance	2"
Turning Radius	35.4"
Overall Length	40"
Overall Width	21"
Overall Height	36"
Seat-to-Floor Height	19" - 21"
Seat-to-Deck Height	16" - 18"
Seat Width	18"
Seat Depth	17"
Backrest Height	15"
Weight of Heaviest pc.	29 lbs.
Total Weight	99 lbs.
Motor	DC 24V, 100W
Controller	PG S-DRIVE 45A
Charger	1.8A off-board
Gradient	6°
Front Wheels	7"
Drive Wheels	8"
HCPCS Code	K0800
Color Options	<div style="display: flex; align-items: center;"> <div style="width: 15px; height: 15px; background-color: red; margin-right: 5px;"></div> Red </div> <div style="display: flex; align-items: center;"> <div style="width: 15px; height: 15px; background-color: blue; margin-right: 5px;"></div> Blue </div>

Features MK batteries



merits



merits Enhancing life through mobility.

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com

