

merits

# ROADSTER S3 S731RS3

Three Wheel Scooter w/Lift Option



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features

- Tiller angle adjustable for ergonomic driving position and easy transferring on and off the unit
- Loop tiller with throttle lever to maneuver by palm or finger
- Transport and store easily with auto latching lockup design
- Detachable battery pack with built in charging port
- Easily disassembles into 5 pieces
- Extra charging port under dashboard for easy plug in
- No tools needed for set-up
- Adjustable seat height and armrest width
- Headlight
- Rear anti-tip wheels prevent a tip-over backwards

## Options

- 40" seat belt
- Multi-accessory receiver adapter
- Oxygen tank holder
- Oxygen tank bag holder
- Crutch and cane holder
- Walker holder
- Medium or Large foldable rear basket
- Cup holder
- Safety flag
- Cover



## Specifications

Product#	S731RS3
Weight Capacity	260 lbs
Drivetrain	RWD
Maximum Speed	4 mph
Per-Charge Range	up to 12 mi.
Battery Size (2 required)	12Ah, ES12 AGM
Ground Clearance	2.2"
Turning Radius	47"
Overall Length	40.5"
Overall Width	25.5"
Overall Height	36.6"
Seat-to-Ground Height	19" - 21"
Seat-to-Deck Height	16" - 18"
Seat Width	18"
Seat Depth	18"
Backrest Height from seat	14.5"
Basket Size	12.5"W x 10"H x 6"D
Weight of Heaviest Piece	29 lbs.
Total Weight	103 lbs.
Motor	DC 24V 180W
Controller	PG S-DRIVE 45A
Charger	1.8A off-board
Gradient	6°
Front Wheels	7"
Drive Wheels	8"
HCPCS Code	n/a
Color Options	<div style="display: flex; align-items: center;"> <div style="width: 15px; height: 15px; background-color: red; margin-right: 5px;"></div> Red         </div> <div style="display: flex; align-items: center;"> <div style="width: 15px; height: 15px; background-color: blue; margin-right: 5px;"></div> Blue         </div>

Features MK batteries



merits



**merits** Enhancing life through mobility.

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0225