

merits

SILVERADO EXTREME S941L

4-Wheel Full Suspension Electric Scooter



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- Travels up to 9.5 mph and can go up to 55 miles on a single charge!
- Sporty streamline design.
- Full suspension system with shock absorbers deliver a smooth and comfortable ride.
- Modern LED lighting system delivers a brilliant light with less power consumption.
- Digital LCD dashboard display: time display, battery gauge, temperature, maintenance messages, speed display, and odometer.
- Extra storage space beneath tiller.
- Automatic safety turning half speed control by sensor.


Accessories

- USB charger
- Oxygen tank holder
- Cup holder & mount

Options

- Therafin positioning products
- Proform headrest with EZ-Clamp bracket

LCD Dashboard Display



LED Headlights & Taillights



Specifications

Product#	S941L
Weight Capacity	450 lbs.
Drivetrain	RWD
Seat Width	20"
Seat Depth	18"
Backrest Height	25"
Seat-to-Floor Height	19"
Seat-to-Ground Height	26.4
Overall Length	61"
Overall Width	29"
Overall Height	57.3"
Turning Radius	82.7"
Ground Clearance	5.1"
Weight of Heaviest Piece	185 lbs.
Total Weight	387 lbs.
Max Speed	9.5 mph
Battery Size (2 required)	100Ah, 12V, Group 27 AGM
Range Per Charge	up to 55 mi.
Motor	DC 24V, 650W
Controller	P&G S-Drive 120A
Charger	8A off-board
Gradient	10°
Wheel Base	45"
Front Wheels	14" Pneumatic Tire
Rear Wheels	16" Pneumatic Tire
Brake	Intelligent, Regenerative, Electromagnetic Brakes.
Color Option	■ Black

Features MK batteries



merits



merits Enhancing life through mobility.

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0325