

merits

# TRAVEL EASE HD P181

Heavy Duty Portable Folding Power Chair

K0824



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

# Features

- Optional weight capacity of 450 lbs.
- 12.5" rubber drive tires for superior traction
- 4.8" ground clearance for easy obstacle clearance
- Foldable for improved storage and transportation

# Specifications

Product#	P181
Weight Capacity	450 lbs., 600 lbs., 700 lbs.
Drive train	RWD
Seat Type	Sling
Seat Width	22"
Seat Depth	18"
Backrest Height	15" - 19"
Seat-to-Floor Height	20"
Overall Length	45"
Overall Width	28"
Overall Height	36"
Turning Radius	25"
Ground Clearance	4.8"
Total Weight	146 lbs. 131 lbs. w/o batteries
Wt of Heaviest Pc.	86 lbs. - 97 lbs.
Max Speed	4 mph
Range Per Charge	up to 18 mi.
Motor	DC 24V / 240W, DC 24V / 950W
Controller	VSI, VR2 70A
Battery Size (2 required)	55Ah, 12V, Group 22NF
Charger	5A off-board
Gradient	10°
Front Wheels	8" foam-filled
Drive Wheels	12.5" foam-filled
Rear Wheels	1.5" anti-tippers
HCPCS Code	K0824
Standard Colors	■ Silvertone

# Accessories

- USB charger
- Oxygen tank holder
- Swing-away joystick
- Cup holder
- Side bag

# Options

- Multiple seat sizes available
- Height adjustable armrests and swing-away footrests
- Elevating leg rests (E0990)



Features MK batteries



merits



**merits** Enhancing life through mobility.

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com

