

merits

TRAVEL EASE P101

Portable Folding Power Chair

K0822



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- Weight capacity of 300lbs.
- 12.5" rubber drive tires for superior traction
- 4.8" ground clearance for easy obstacle clearance
- Foldable for improved storage and transportation

Accessories

- USB charger
- Oxygen tank holder
- Swing-away joystick
- Cup holder
- Side bag

Options

- Height adjustable armrests and swing-away footrests
- Elevating leg rests (E0990)



Specifications

Product#	P101
Weight Capacity	300 lbs.
Drive train	RWD
Seat Type	Sling
Seat Width	18" - 20"
Seat Depth	16"
Backrest Height	15" - 19"
Seat-to-Floor Height	20"
Overall Length	43"
Overall Width	24" - 26"
Overall Height	36"
Turning Radius	23"
Ground Clearance	4.5"
Total Weight	150 lbs. - 155 lbs.
Wt of Heaviest Pc.	86 lbs. - 97 lbs.
Max Speed	4 mph
Range Per Charge	up to 18 mi.
Motor	DC 24V / 180W, DC 24V / 630W
Controller	PG VSI 50
Battery Size (2 required)	33Ah, U1 AGM
Charger	5A off-board
Gradient	10°
Front Wheels	8" foam-filled
Drive Wheels	12.5" foam-filled
Rear Wheels	1.5" anti-tippers
HCPCS Code	K0822
Standard Colors	■ Silvertone

Features MK batteries



merits



merits Enhancing life through mobility.

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0125