



# ULTRA HD P331

Heavy Duty MWD Powerbase

K0837



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features

- Mid-wheel drive design provides outstanding maneuverability
- Adjustable seat width 20"-24" and seat depth 16"-22" for custom fit
- Optimized seating stability with 4-post seat design
- Tool-free battery box provides easy maintenance
- Front and rear suspension offers all day comfort and easy obstacle clearance
- Fixed back angle: 90°, 95°, 100°, 105° and 110°
- Dual in-line motors for enhanced efficiency, torque, range, and performance
- Unoccupied transit tie down brackets included

## Accessories

- USB charger
- Oxygen tank holder bag

## Options

- Tilt 0° - 50°
- Recline 90° - 110° fixed
- Articulating Footplate
- VR2 electronics
- Swing away mounting brackets and mounts
- Multiple Joystick options
- Power Seat Tilt/ Angle-Adjustable Back
- Seat to Floor Height options
- Backrest Width 12" -24"
- Positioning Components and Headrests

## Specifications

Product#	P331
Weight Capacity	450 lbs
Drivetrain	MWD
Maximum Speed	5 mph
Per-Charge Range	up to 8 miles / 19 miles
Battery Size (2 required)	55Ah, 22NF AGM 51Ah, 22NF GEL 33Ah, U1 AGM
Ground Clearance	2.875"
Turning Radius	21"
Overall Length	35" / 45" with footplate
Overall Width	27"
Overall Height	42"
Seat-to-Ground Height	17.5" - 18.5"
Seat Width	20" - 24"
Seat Depth	16" - 22"
Backrest Height	21" - 23"
Total Weight	288lbs. - 310 lbs
Motor	2-Pole/200W
Controller	PG VR2-A2 60A
Charger	6 Amp
Gradient	6°
Front Wheels	6" x 2"
Drive Wheels	12" x 3"
Rear Wheels	6" x 2"
HPCSCCode	K0837
Color Options	<input type="checkbox"/> White <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> Charcoal Gray



Features MK batteries



 **merits** *Enhancing life through mobility.*

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0225