



# ULTRA 3A P325

Tilt | Elevate | Articulating Foot Platform

K0835

**TIME Business**  
**NEW PRODUCT AWARDS**

★ **2024** ★



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features

- Center of gravity tilt 0° - 50° posterior
- Articulating Foot Platform (AFP) with simultaneous extension and articulation
- Tension-adjustable front and rear suspension
- Adjustable seat width and depth from 16" - 22"; depth adjustable in 1" increments
- 24" base width and seat-to-floor height as low as 16.5"
- Rated up to 19 miles per charge

## Options

- 22NF battery upgrade
- 12" elevating seat
- Power articulating foot platform
- Elevating leg rests or swing-away footrests
- Swing-away joystick
- Oxygen tank holder
- USB charger



Group 2  
Single Power  
3 Actuators

## Specifications

Product#	P325
Weight Capacity	300 lbs.
Drivetrain	MWD
Maximum Speed	5 mph
Per-Charge Range	up to 19 mi.
Battery Size (2 required)	55Ah, 22NF AGM 51Ah, 22NF GEL 33Ah, U1 AGM
Ground Clearance	2.5"
Turning Radius	26"
Overall Length	36"
Overall Width	24"
Suspension	Front and Rear
Seat-to-Floor Height	Power tilt or pan: 16.5"-19.5", Captain: 20"-23"
Seat Width	16"-22", adjusted in 1" increments
Seat Depth	16"-22", adjusted in 1" increments
Backrest Height	18"-28"
Total Weight	315 - 384 lbs.
Motor	2-Pole/240W
Controller	P&G VR2-A2 60A
Charger	6 Amp
Gradient	6°
Front Wheels	Anti-Tipper
Drive Wheels	14" x 3" foam-filled
Casters	6" x 2" foam-filled
HCPCSCode	K0835
Color Options	<input type="checkbox"/> White <input checked="" type="checkbox"/> Red <input checked="" type="checkbox"/> Blue <input checked="" type="checkbox"/> Charcoal Gray

Features MK batteries



**merits** Enhancing life through mobility.

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com

