



# VISION ULTRA X P325X

MWD Compact Single Power Base

K0837



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

# Features

- 330 lb. capacity
- 14" drive tires
- Tension-adjustable front and rear suspension
- Unoccupied tie down restraints
- Multiple seating options, power tilt, rehab seat frame, captain or solid seat pan
- Multiple seat-to-floor height options: 16.5" - 19.5"
- Width and depth-adjustable seat; 16" - 22" in 1" increments
- Red, white, blue and charcoal gray color options

# Accessories

- USB charger
- Weather cover

# Options

- 22NF battery upgrade
- Elevating leg rests
- Swing-away foot-rests
- Swing-away joystick
- Power articulating foot platform



# Specifications

Product#	P325X
Weight Capacity	330 lbs.
Drivetrain	MWD
Maximum Speed	5 mph
Per-Charge Range	8.38 mi.
Battery Size (2 required)	55Ah, 22NF AGM 51Ah, 22NF GEL 33Ah, U1 AGM
Ground Clearance	2.5"
Turning Radius	27"
Overall Length	46.6"
Overall Width	25.75
Suspension	Front and Rear
Seat-to-Floor Height	Power tilt 16.5" - 19.5"
Seat Width	16" - 22"
Seat Depth	16" - 22"
Backrest Height	21" - 23"
Total Weight	355.5 lbs.
Motor	2-Pole/240W
Controller	PG VR2-A2 60A
Charger	6 Amp
Gradient	6°
Front Wheels	Anti-Tipper
Drive Wheels	14" x 3" foam-filled
Casters	6" x 2" foam-filled
HCPCS Code	K0837
Color Options	<div style="display: flex; flex-direction: column; gap: 5px;"> <div><span style="display: inline-block; width: 15px; height: 15px; background-color: red; border: 1px solid black;"></span> Red</div> <div><span style="display: inline-block; width: 15px; height: 15px; background-color: white; border: 1px solid black;"></span> White</div> <div><span style="display: inline-block; width: 15px; height: 15px; background-color: blue; border: 1px solid black;"></span> Blue</div> <div><span style="display: inline-block; width: 15px; height: 15px; background-color: gray; border: 1px solid black;"></span> Charcoal Gray</div> </div>

Features MK batteries



**merits** *Enhancing life through mobility.*

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com

