



AXCEL NARROW P335N

FWD Powerbase

K0856 / K0861



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- 22.5" wide front wheel drive chair is easy to maneuver in tight spots, narrow hallways and door frames
- Tight 90° turns
- Navigating over uneven terrain is easy with FWD 14" drive tires, maintaining contact.
- Independent hydraulic, high-performance suspension
- Large 9" rear casters to navigate obstacles and rough terrain
- Bright LED lights, included
- Tool-free, rear battery access
- Power recline adjustments, 85° - 160°
- Powered by MK batteries, rated up to 19 miles per charge

Accessories

(Based on chair configuration)

- USB charger
- Oxygen tank holder
- Cup holder & mount

Options

- ASL alternative drive
- Therafin positioning products
- Power articulating footplate
- Push handles
- Proform headrest with EZ-Clamp bracket
- 12" seat elevation

Specifications

Product#	P335N
Weight Capacity	300 lbs.
Drivetrain	FWD
Seat Width	14 - 20"
Seat Depth	14 - 20"
Backrest Height	18 - 22"
Seat-to-Floor Height	17.5 - 19.5"
Seat Elevation (Optional)	12"
Overall Length	39.6"
Overall Width	22.5"
Overall Height	41.5"
Turning Radius	26" w/o footplate
Ground Clearance	3"
Total Weight	350-385 lbs.
Max Speed	5.2 mph
Battery Size (2 required)	50Ah, M50 AGM
Range Per Charge	up to 19 mi.
Motor	2-Pole/222W
Controller	R-NET 90A
Charger	6A off board
Gradient	7.5°
Front Wheels	Anti-Tipper
Drive Wheels	14"x 3" hybrid gel
Rear Wheels	9" hybrid gel
HCPCode	K0856/K0861
Color Options	<div> <div>Black Noir</div> <div>Dove Gray</div> <div>Sapphire Blue</div> <div>Ruby Red</div> <div>Camo Green</div> </div>



Features MK batteries



merits Enhancing life through mobility.

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0425