


GEL


The PC Gel Cushion features an ergonomic contour with a gel insert overlay, offering an excellent entry-level option for new Precision Comfort users. It offloads coccyx/sacrum and supports healthy pelvic positioning for lower extremity stability. The gel overlay enhances immersion and support, ensuring comfort during prolonged sitting.




Key features / Clinical justification



Improves pelvic and lower extremity balance for postural stability



Rear cut-out for effective coccyx/sacrum offloading to prevent pressure injury



Immersive gel overlay for effective weight redistribution further reducing pressure injury risks

HPCPS E2603, E2604

Custom sizing available upon request

Specification

450lb Max Load

Sizes	14"W	16"W	18"W	20"W	22"W	24"W
14"D	14x14	16x14	-	-	-	-
16"D	14x16	16x16	18x16	20x16	-	-
18"D	14x18	16x18	18x18	20x18	22x18	24x18
20"D	-	16x20	18x20	20x20	22x20	24x20
22"D	-	-	-	-	22x22	24x22
24"D	-	-	-	-	-	24x24