

# HEADREST


The PC headrest combines adaptive support with universal wheelbase compatibility. The multi-layer foam pad balances sustained comfort with robust support, while precision positioning is achieved through a stainless-steel gear system. Available in mid-line or multi-axis configurations, the 6" stem provides extensive adjustment range.

The EZClamp™ bracket features single-latch, tool-free removal for simplified handling. This system delivers optimal positioning with enhanced usability for diverse clinical needs.




EZClamp™  
bracket

## Key features / Clinical justification



Layered padding enhances comfort and protects skin integrity



Prevents postural abnormalities progression that may affect digestion and respiration

Reduces effort in maintaining an upright head position, supporting engagement with the environment

HPCPS E0955, E1028

## Specification

- 6" Hardware stem
- 2 Removable hardware options: mid-line or multi-axis
- Toolless EZClamp™ bracket