

LATERAL TRUNK SUPPORT

The PC lateral trunk support accommodates both low and high muscle tone with engineered durability for clients with spasticity. Its linked rotation mechanism ensures optimal pad contact by conforming to individual body shapes, while enabling precise adjustments for postural asymmetries.




Quick release
button


Swing away for easy side transfer.




Key features / Clinical justification



Correct or support body asymmetries, preventing further spinal deformities



Enhances core stability to improve upper motor function



Upper trunk support improves respiration and digestion

HPCPS E1028 / E1034 for lateral hardware, E0956 for positioning pads (Coming Soon)

Specification

- 1.25" extension length
- ~1.25" offset
- Withstands user weight up to 300lb weight capacity

Available Pad sizes

3"x5"
4"x6"
4"x8"
4x12"