


# LITE

The PC Lite is a high-performance general use cushion for users who can independently perform weight shifts. Its ergonomic contour enhances lower extremity stability and reduces fatigue. A rear concave channel offloads coccyx pressure, while the moisture-wicking cover keeps the user dry. Its lightweight design allows easy transfers between chairs.




## Key features / Clinical justification



Comfort and support  
beyond sling upholstery

Ergonomic design for  
lower extremity stability



Coccyx offloaded with cut  
out design at the rear of  
cushion

HPCPS E2601, E2602

### Specification

Sizes	16"W	18"W	20"W	22"W	24"W
16"D	16x16	18x16	20x16	-	-
18"D	16x18	18x18	20x18	22x18	24x18
20"D	16x20	18x20	20x20	22x20	24x20
22"D	-	-	-	22x22	24x22
24"D	-	-	-	-	24x24

450lb Max Load