



A DIVISION OF MERITS  
HEALTH PRODUCTS

FOR IMMEDIATE RELEASE

1/3/2025

## PILOT UNVEILS REVOLUTIONARY SPACE-SAVING STAIRLIFT: THE PILOT FLIGHT

Fort Myers, FL – 1/3/2025 – Pilot Stairlifts, a global leader in advanced leader mobility solutions, is excited to announce the debut of the Pilot Flight, an innovative stairlift that redefines convenience, safety, and space efficiency for individuals with limited mobility. Designed specifically for straight staircases, the Pilot Flight combines cutting-edge technology with user-centric features to deliver unparalleled performance and adaptability.

The Pilot Flight stairlift sets itself apart with its innovative enclosed gear rack and slim carriage design, which allows the lift to fold compactly against the wall when not in use. This ingenious feature not only conserves valuable living space but also maintains a clean and open atmosphere within the home. Users can experience a seamless transition as the stairlift unfolds and descends from its compacted position, offering reliable support up and down the stairs. The Pilot Flight stairlift is a testament to Pilot Stairlifts' dedication to merging state-of-the-art technology with practical design, ensuring a product that meets and exceeds the needs of individuals with mobility challenges.

"We are excited to introduce the Pilot Flight, which provides independence and freedom of movement for those with mobility challenges," said Andrew Maskery President of Merits. "Our goal is to improve lives through high-quality mobility solutions, and the Pilot Flight stairlift does just that in an elegant, space-efficient design."

### Pilot Flight Key Features:

- Space-saving design: enclosed track spacing measures only 4.6" from the wall.
- Functionality: equipped with a diagnostic display and easy access panels for fast, accurate servicing.
- Smooth operation: soft start and stop mechanisms provide a gentle and comfortable ride for users.
- Enhanced comfort: features a swivel seat with a locking catch, facilitating easy and secure entry and exit.
- Easy to maintain: diagnostic display and quick access panels.
- Safety first: the stairlift includes a comprehensive suite of safety features such as a five-directional safety sensor, armrest position sensors, a seat position sensor, an over-speed governor, an electromechanical brake, a seatbelt, and a key switch.
- User capacity: accommodates individuals weighing up to [300] lbs.
- Flexible design: choose installation of the track for the left/right side of the staircase.

### Merits Health Products:

*Pilot Stairlifts, a division of Merits Health Products specializes in giving people with mobility issues the freedom and independence to remain in the home they love by providing stylish, elegant, and affordable solutions to getting up and down the stairs effortlessly.*

*Merits USA, established in 1995, is a subsidiary of Merits Health Products headquartered in Fort Myers, Florida. With a mission to enable independence and accessibility for all, Merits provides quality mobility solutions designed to support people in fully experiencing all that life has to offer. Merits maintains a global presence with distribution centers located on both U.S. coasts, allowing them to effectively serve client needs. Manufactured in the United States, Merits Health Products is committed to improving lives worldwide daily through the development of high-quality power wheelchairs, scooters, and stairlifts.*

### Media Contact:

Olesya Timofeeva

Marketing Coordinator

(239) 427-4088 EXT: 0368

[olesya.timofeeva@meritsusa.com](mailto:olesya.timofeeva@meritsusa.com)

[www.meritsusa.com](http://www.meritsusa.com)

