

merits

REGAL P310

Compact RWD Powerbase

K0822 / K0823



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- 10" knobby drive tires for superior traction
- 3" ground clearance for easy obstacle clearance
- Dynamic LiNX controller for easy programming and diagnostics
- 16", 18" or 20" ergonomic, semi-reclining seats: up to 135°

Accessories

- USB charger
- Oxygen tank holder
- Walker holder
- Cover
- Cup holder
- Mounting bracket for cup holder

Options

- Medium or Large foldable rear basket
- Side bag
- Large footplate
- Elevating legrests
- Captain or Pan seat

Specifications

Product#	P310
Weight Capacity	300 lbs.
Drivetrain	RWD
Seat Type	Captain or Pan
Seat Width	16", 18" or 20"
Seat Depth	16" or 18"
Backrest Height	18" or 20"
Seat-to-Floor Height	20"-22" (captain), 16.5"-18.5" (pan)
Overall Length	39"
Overall Width	24"
Overall Height	18" or 20"
Turning Radius	21"
Ground Clearance	3"
Total Weight	197.5 lbs.
Max Speed	5 mph
Range Per Charge	up to 18 mi.
Motors (x2) Left & Right	M13CM-4 LH M13CM-4 RH
Controller	Controller Dynamic LiNX
Battery Size (2 required)	33Ah, U1 AGM
Charger	5A off-board
Gradient	12°
Front Wheels	8" foam filled tire
Drive Wheels	10" foam filled tires
Rear Wheels	3" anti-tippers
HCPSC Codes	K0822 / K0823
Color Options	<div><div></div> Red</div> <div><div></div> Blue</div>



Features MK batteries




merits



merits Enhancing life through mobility.

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0425