



# ULTRA X P325X

MWD Compact Single Power Base

K0837



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features

- 330 lb. capacity
- 14" drive tires
- Tension-adjustable front and rear suspension
- Unoccupied tie down restraints
- Multiple seating options, power tilt, rehab seat frame, captain or solid seat pan
- Multiple seat-to-floor height options: 16.5" - 19.5"
- Width and depth-adjustable seat; 16" - 22" in 1" increments
- Red, white, blue and charcoal gray color options

## Accessories

- USB charger
- Weather cover

## Options

- 22NF battery upgrade
- Elevating leg rests
- Swing-away foot-rests
- Swing-away joystick
- Power articulating foot platform



## Specifications

Product#	P325X
Weight Capacity	330 lbs.
Drivetrain	MWD
Maximum Speed	5 mph
Per-Charge Range	12 mi. / 19 mi.
Battery Size (2 required)	34Ah/U1 or 55Ah/22NF
Ground Clearance	2.5"
Turning Radius	26"
Overall Length	36"
Overall Width	24"
Suspension	Front and Rear
Seat-to-Floor Height	Power tilt 16.5" - 19.5"
Seat Width	16" - 22"
Seat Depth	16" - 22"
Backrest Height	21" - 23"
Total Weight	315 lbs. - 385 lbs.
Motor	2-Pole/240W
Controller	PG VR2-A2 60A
Charger	6 A off board
Gradient	6°
Front Wheels	6"x 2" solid
Drive Wheels	14"x 3" foam-filled
Casters	6"x 2" solid
HCPCS Code	K0837
Color Options	<div> <div style="display: inline-block; width: 15px; height: 15px; background-color: red; margin-right: 5px;"></div> Red                 <div style="display: inline-block; width: 15px; height: 15px; background-color: white; border: 1px solid black; margin-right: 5px;"></div> White                 <div style="display: inline-block; width: 15px; height: 15px; background-color: blue; margin-right: 5px;"></div> Blue                 <div style="display: inline-block; width: 15px; height: 15px; background-color: gray; margin-right: 5px;"></div> Charcoal Gray             </div>

Features MK batteries



**merits** Enhancing life through mobility.

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0425