



ULTRA 3A P325

Tilt | Elevate | Articulating Foot Platform

K0835

HME Business
NEW PRODUCT
AWARDS

★ **2024** ★



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- Center of gravity tilt 0° - 50° posterior
- Articulating Foot Platform (AFP) with simultaneous extension and articulation
- Recline angle 90° - 110°
- Tension-adjustable front and rear suspension
- Adjustable seat width and depth from 16" - 22"; depth adjustable in 1" increments
- 24" base width and seat-to-floor height as low as 16.5"
- Rated up to 19 miles per charge with 22NF batteries

Options

- 22NF battery upgrade
- 12" elevating seat
- Power articulating foot platform
- Elevating leg rests or swing-away footrests
- Swing-away joystick
- Oxygen tank holder
- USB charger



Group 2
Single Power
3 Actuators

Specifications

Product#	P325
Weight Capacity	300 lbs.
Drivetrain	MWD
Maximum Speed	5 mph
Per-Charge Range	up to 19 mi.
Battery Size (2 required)	55Ah, 22NF AGM 51Ah, 22NF GEL 33Ah, U1 AGM
Ground Clearance	2.5"
Turning Radius	26"
Overall Length	36"
Overall Width	24"
Overall Height	40"
Elevation	12"
Suspension	Front and Rear
Seat-to-Floor Height	Power tilt or pan: 16.5"-19.5", Captain: 20"-23"
Seat Width	16"-22", adjusted in 1" increments
Seat Depth	16"-22", adjusted in 1" increments
Backrest Height	21"-23"
Total Weight	384 lbs.
Motor	DC 24V, 240V
Controller	P&G VR2-A2 60A
Charger	6A off board
Gradient	6°
Front Wheels	6" x 2" foam-filled
Drive Wheels	14" x 3" foam-filled
Casters	6" x 2" foam-filled
HCPSCCode	K0835
Color Options	<input type="checkbox"/> White <input checked="" type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> Charcoal Gray

Features MK batteries



merits Enhancing life through mobility.

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0425