

merits

VISION SUPER P327

Super Heavy Duty MWD Powerbase

K0824 / K0825 / K0830 / K0831



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- True mid-wheel: 21" turning radius
- Robust independent suspension can absorb up to 350 lbs. on the footplate without tipping forward (does not apply to lift option)
- Dynamic LiNX & Shark controller for convenient Bluetooth programming and diagnostics
- Premium, semi-reclining ergonomic seat up to 135°
- 18" - 22" customizable seat depth

Accessories



- Unoccupied transit tie down kit
- Cup holder
- USB charger
- Oxygen tank holder
- Walker holder
- Foldable rear basket
- Safety flag
- Side under, arm bag

Options

- Optional, 20" - 24" seat width
- 10" scissor power seat lift
- Elevating leg rests (K0195)
- Captain or Pan powered seat option



Specifications

Product#	P327/P3274 w/lift
Weight Capacity	450 lbs., 400lbs. w/lift
Drivetrain	MWD
Seat Type	Captain or Pan
Seat Width	20" - 24"
Seat Depth	18" - 22"
Backrest Height	20" - 22"
Seat-to-Floor Height	22" 23"/24" - 34" w lift
Seat Elevation	10" - P3274 w/lift
Overall Length	35.4"
Overall Width	27.2"
Overall Height	48"
Turning Radius	21"
Ground Clearance	2.2"
Total Weight	243 lbs., 273 lbs.
Wt of Heaviest Pc.	121 lbs., 151 lbs.
Max Speed	5 mph
Range Per Charge	up to 20 mi.
Motor	DC 24V / 300W x 2pcs
Controller	Dynamic Shark 50A / PG VR2 60A (Lift)
Battery Size (2 required)	55Ah, 22NF AGM
Charger	5A off-board
Gradient	8°
Front Wheels	6" - solid
Drive Wheels	12" - flat-free
Rear Wheels	6" - solid
HCPSC Codes	K0824 / K0825 K0830 / K0831
Color Options	 Red  Black

Features MK batteries



merits



 **Enhancing life through mobility.**

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0425