

merits

GEMINI P301

Heavy Duty RWD Powerbase

K0824 / K0825



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- Weight capacity of 450 lbs.
- Ergonomic semi-reclining seat that reclines up to 135°
- 14" rubber drive tires for superior traction
- True RWD and 35" turning radius
- Optional 5" power elevating seat
- Rated up to 5 mph and range up to 20 miles

Accessories



- Oxygen tank holder
- Cane & crutch holder
- Walker holder
- Foldable rear basket
- Cup holder
- Side bag
- Cover
- USB charger

Options

- Captain or Pan powered seat option
- 20" and 22" seat width option
- Elevating leg rest available (K0195)
- 5" power lift option, available
- Swing-away joystick
- XL footplate



Specifications

Product#	P301
Weight Capacity	450 lbs.
Drivetrain	RWD
Seat Type	Captain or Pan
Seat Width	21.75"
Seat Depth	18" - 22"
Backrest Height	21.5"
Seat-to-Floor Height	23"
Overall Length	29.5"
Overall Width	28"
Overall Height	42.75"
Turning Radius	35"
Ground Clearance	2.5"
Total Weight	242 lbs.
Wt of Heaviest Pc.	126 lbs.
Max Speed	5 mph
Range Per Charge	up to 20 mi.
Motor	DC 24V / 240W
Controller	P&G VR2 60A
Battery Size (2 required)	50Ah, 22NF, 12V
Charger	5A off-board
Gradient	10°
Front Wheels	8" - foam-filled
Drive Wheels	14" - foam-filled
Rear Wheels	4" anti-tippers
HCPCS Codes	K0824 / K0825
Color Options	 Candy Apple Red  Metallic Silver Blue

Features MK batteries



 **Enhancing life through mobility.**

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0525