

ROADSTER 4 S741A

Four Wheel Scooter

K0800



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- 4-wheel design provides increased stability
- Easily disassembles into 5 pieces for transport and storage.
- 18" wide deluxe stadium seat
- Detachable battery pack
- Tiller angle adjustable for ergonomic driving position and easy transferring on and off the unit
- Loop tiller with throttle lever to maneuver by palm or finger
- Transport and store easily with auto latching lockup design
- Adjustable seat height and armrest width
- No tools needed for set-up

Options

- 40" seat belt
- Multi-accessory receiver adapter
- Oxygen tank holder
- Oxygen tank bag holder
- Crutch and cane holder
- Walker holder
- Medium or Large foldable rear basket
- Cup holder
- Safety flag
- Cover











Specifications

S741A

300 lbs

RWD

4 mph

43.3"

30"

21"

36"

18" 17"

15"

6°

8" 8"

K0800

Red

36 lbs.

106 lbs.

DC 24V, 100W

PG S-DRIVE 45A

1.8A off-board

19" - 21"

16" - 18"

up to 6 miles

12Ah, ES 12 AGM

Product#

Drivetrain

Weight Capacity

Maximum Speed

Per-Charge Range

Ground Clearance

Turning Radius

Overall Length

Overall Width

Overall Height

Seat Width

Seat Depth

Backrest Height

Total Weight

Motor

Controller

Charger

Gradient

Front Wheels

Drive Wheels **HCPCS** Code

Color Option

Seat-to-Floor Height

Seat-to-Deck Height

Weight of Heaviest pc.

Battery Size (2 required)















Enhancing life through mobility.









