

merits

VISION CF P322

FWD Powerbase

K0815 / K0816 / K0820 / K0821 / K0822 / K0823



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- A true FWD with a 17" turn radius
- 9" knobby drive tires for superior traction
- Dynamic LiNX controller for convenient Bluetooth programming and diagnostics
- Premium, semi-reclining ergonomic seat up to 135°
- 16" - 20" customizable seat depth in 1" increments

Accessories

- Cup holder
- USB charger
- Oxygen tank holder
- Walker holder
- Foldable rear basket
- Safety flag
- Side under, arm bag

Options

- Optional, 16" - 20" seat width
- Elevating leg rests (K0195)
- Captain or Pan powered seat option



Specifications

Product#	P322
Weight Capacity	300 lbs.
Drive train	FWD
Seat Type	Captain or Pan
Seat Width	16" - 20"
Seat Depth	16" - 20"
Backrest Height	18" - 20"
Seat-to-Floor Height	16.5" - 22"
Overall Length	33"
Overall Width	24"
Overall Height	37"
Turning Radius	17"
Ground Clearance	2.4"
Total Weight	154 lbs.
Wt of Heaviest Pc.	62.7 lbs.
Max Speed	4 mph
Range Per Charge	up to 18 mi.
Motor	DC24V / 130W
Controller	Dynamic LiNX 40A
Battery Size (2 required)	33Ah, U1 AGM
Charger	4A off-board
Gradient	6°
Front Wheels	3" anti-tippers
Drive Wheels	9" foam-filled
Rear Wheels	6" foam-filled
HCPCS Code	K0815 / K0816 K0820 / K0821 K0822 / K0823
Color Option	 Candy Apple Red

Features MK batteries



merits



 **Enhancing life through mobility.**

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0425