

SIZE MATTERS!

AXCEL NARROW

Avid Rehab



14"
FWD DRIVE



22.4"
FWD BASE

6.45 MPH
MAX. SPEED

300 LBS
WEIGHT CAPACITY

12" LIFT
ELEVATION

www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



©Merits Health Products, Inc. | 4245 Evans Avenue Fort Myers, FL 33901

AXCEL NARROW P335N

With **FREQUENCY** Seating

K0856 / K0861



Specifications

Product name	Axcel Narrow
Model #	P335N
Group	Group 3
HCPCS Code	K0856 / K0861
Weight Capacity	300 lbs.
Drivetrain	FWD
Seat Type	FREQUENCY Seating
Seat Width	14 - 20"
Seat Depth	14 - 20"
Seat-to-Floor Height	17.5 - 19.5"
Backrest Adjustment Range	23 - 30"
Elevating Seat	12"
Recline Angles	85° - 160°
Motor	2-Pole / 222W
Max Speed*	6.45 mph
Controller	R-NET 90A
Turning Radius	26" w/o footplate
Ground Clearance	3"
Gradient	7.5°
Battery Size (2)	50Ah, M50 AGM
Per Charge Range	up to 12 mi.
Charger	8A off board
Front Wheels	Anti-Tipper
Drive Wheels	14"x 3" hybrid gel
Rear Wheels	9" hybrid gel
Overall Length	39.6"
Base Width	22.4"
Overall Height	41.5"
Total Weight*	350-430 lbs.
Standard Colors	Black Gray Matte Blue Red Matte Green

*Depends on configuration

Features

- Industry's first 22.4-inch front-wheel-drive power wheelchair, engineered for unmatched maneuverability without sacrificing performance
- Up to 6.45 mph* maximum speed allows for increased mobility and independence
*With Advanced Stability Control
- Independent hydraulic suspension and 14-inch front drive wheels ensure a smooth, stable ride, while large 9-inch rear casters handle obstacles with ease.
- FREQUENCY Seating allows for adjustable power recline with 6" shear reduction for a comfortable reclining experience
- Tool-free battery access and modular shroud for easy maintenance
- Bright LED front and rear lights included for safety

Accessories

(Based on chair configuration)

- USB charger
- Oxygen tank holder
- Cup holder & mount

Options

- ASL alternative drive
- Precision Comfort positioning products
- Power articulating footplate
- Push handles
- ProForm headrest with EZ-Clamp bracket
- 12" seat elevation



merits Enhancing life through mobility.

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com

Features MK batteries



REV.0926